



MEDICAL DEVICES AND P.P.E. PRODUCTION

COMPANY & PRODUCTION

iVision Health is an industry specialized in manufacturing of medical devices and personal protection equipment (P.P.E.)

Our production plants are located in Italy, in Martignacco (UD) and in San Dorligo della Valle (TS), in a 12,000 square meters space.

Our strengths:

- Production capacity of 15,000,000 surgical masks per month and more than 1.000,000 filtering facial-half-masks (FFP) per month
- Latest generation machinery
- Use of the **highest quality Italian fabrics**
- Process and product quality certifications and authorizations: ISO, CE, ISS
- A complete range of products in response to market needs




MADE IN ITALY PRODUCTION

iVision Health uses latest **CE certificated** machinery generation, characterized by high precision and speed.

Our production process protocol provides stringent hygiene and cleaning standards, sanification of environment, plants and production lines, in order to guarantee the quality of our products.

The **iVision Health** factory consists of:

- Dedicated lines for mask production
- Packaging automation for different formats
- Trained and qualified internal staff



ITALIAN FABRICS


iVision Health products are made with Italian materials of the highest quality. We carry out careful and continuous research in this regard

Our pluses are:

- Very high bacterial filtration capacity
- Maximum comfort and softness for extended daily usage
- Soft fabrics free of irritants
- Constant security assured by in-depth periodic laboratory test

The raw materials we use are as follows:

- Soft spunbond non-woven fabric
- Particulate and antibacterial filter (meltblown / SMS)
- Anti-bacterial, anti-static, breathable, ecological Air Cotton fabric



CERTIFICATIONS



iVision Health is not only a certified company but it also holds all the necessary authorizations in order to produce surgical masks, according to the standards required by law:

- ISO 9001:2015 (quality certification) for the design, the production and the marketing of surgical masks
- ISO 13485:2016 (medical device quality certification) for the design, production and marketing of surgical masks
- ISS Authorization Istituto Superiore di Sanità -COV1297 and COV1303 reference www.iss.it
- CE Certification, according to 93/42 CEE directive www.salute.gov.it









Questo certificato è stato rilasciato a

IVISION S.R.L.

Sede legale: VIA DEI SOLTERI, 105 - 38121 TRENTO, Italy Sede coprativa: VIA J.RESSEL, 2 - 34018 SAN DORLIGO DELLA VALLE DOLINA (TS)

come riconoscimento del Sistema Gestione Qualità aziendale in conformità alle Norme

ISO 13485-2016

Le attività coperte da questo certificato sono

Progettazione, Produzione e Commercializzazione di Mascherine Chirurgiche di Tipo II

Numero Del Certificato:

Data d'emissione (original

Data d emissione

Emissione No:

Data di Scadenza:

15 Ottobre 2023

nesso da: Ommud

A nome dal Saheme Manager









SURGICAL MASK

Surgical masks are intended to prevent the contamination of the environment by the wearer, as they limit the transmission of infectious agents. According to Legislative Decree 24 February 1997, Hen.46 and subsequent amendments. (source: http://www.salute.gov.it), surgical masks are included among medical devices category.

iVision Health produces surgical masks of different types, in accordance with the best technical characteristics of the tissues.

iVision Health masks guarantee:

- Liquid splash resistance
- Breathability
- Bacterial filtration efficiency
- Cleaning by microbes



FILTERING FACIAL HALF MASKS

The **filtering facial half masks** (FFP2 and FFP3) are used in hospital and care settings to protect the user from external agents (including the transmission of infections by droplets and aerosols), are certified in accordance with the provisions of Legislative Decree n. 475/1992. They are **Individual Protection Devices** (source: http://www.salute.gov.it).

iVision Health produces different types of filtering face half masks, guaranteeing respiratory protection in accordance to UNI EN 149:2009 technical standard:

- **Dust filter**
- **Guarantee on efficiency characteristics**
- **Breathability**
- **Stability of the structure**



PRODUCTS





SURGICAL MASK









iVi98 SURGICAL MASK

iVi98S SURGICAL MASK

iVilight / iViWhite SURGICAL MASK

Type II high filtering mask

Outer Layer: 25 gr. Filter Layer: 22 gr. Outer Layer: 25 gr.

B.F.E. ≥ 98%

Type II high filtering mask

Outer Layer: 25 gr. Filter Layer: 35 gr. Outer Layer: 25 gr.

B.F.E. ≥ 98%

Type II high filtering mask

Outer Layer: 25 gr. Filter Layer: 35 gr.

B.F.E. ≥ 98%







PRODUCTS





SURGICAL MASK









iViBlue SURGICAL MASK



iViColor SURGICAL MASK

Type II high filtering mask

Outer Layer: 25 gr. Filter Layer: 35 gr. Outer Layer: 25 gr.

B.F.E. ≥ 98%



Outer Layer: 20 gr. Filter Layer: 25 gr. Outer Layer: 25 gr.

B.F.E. ≥ 98%

Type II high filtering mask

Outer Layer: 25 gr. Filter Layer: 35 gr. Outer Layer: 25 gr.

B.F.E. ≥ 98%







PRODUCTS





FILTERING MASKS





FFP2

Filtering Facial Half Masks
Outer Layer: 50 gr. (SPUN)
Filter Layer: 25 gr. (MELT)
Filter Layer: 25 gr. (MELT
Inner Layer: 40 gr. (AIR COTTON)

Outer Layer: 30 gr. (SPUN)

B.F.E. ≥ 95%



FFP3

WORK IN PROGRESS

B.F.E. ≥ 98%

AVAILABLE IN APRIL 2021



MODELS WITH VALVE


WORK IN PROGRESS

B.F.E. ≥ 98%

AVAILABLE IN APRIL 2021

COMMERCIAL MATRIX





SOME OF OUR CLIENTS





















































MASK GENERAL LEGEND





SURGICAL MASKS

The main intended use of **surgical masks** is to **protect people around us** from infectious agents. They can also be worn by patients and other people in order to reduce the risk of spreading infections, particularly in epidemic or pandemic situations.

Surgical masks are classified worldwide as Type I and II:

- Type I masks have a filtration capacity of (B.F.E) ≥ 95%
- Type II masks have a filtration capacity of (B.F.E) ≥ 98%



FILTERING MASKS (FFP1,FFP2,FFP3)

The main function of **Filtering Masks** with valve is to **protect only the wearer** as they filter inwards but not outwards, therefore, an infected person can infect the others. Filter masks are mainly used in professional field (doctors in infectious wards, etc.). They are divided into **three classes of FFP** protection (the acronym stands for filtering face piece, in Italian "particle filtering facepiece")

- FFP1 have a filtration capacity of at least 80%
- FFP2 have a filtrarion capacity of at least 95%
- FFP3 have a filtration capacity of at least 99%





IVISION HEALTH by IVISION GROUP

Office

Via Spilimbergo 154, Martignacco (UD) Tel +39 0432 1483803

Production Plant UD

Via Spilimbergo 154, Martignacco (UD)

Production Plant TS

Via J. Ressel 2, San Dorligo della Valle (TS)

info@ivisionhealth.it www.ivisionhealth.it www.ivision-group.it